


# HelXAS: Design and function

Ari-Pekka Honkanen


# GOAL: THE ABSORPTION COEFFICIENT



Attenuation of photons connected to the properties of matter

*How to measure in home lab?*

# Parts needed:




1. Source of X-rays
2. Device for wavelength separation
3. Photon detector

Image source: [www.ikea.com](http://www.ikea.com)

+ some control and acquisition electronics..



# HeIXAS



# X-RAY TUBE - POLYCHROMATIC SOURCE






Image source: <https://commons.wikimedia.org/wiki/File:XrtSpectrum.jpg>


# Ag TUBE SPECTRUM




- HelXAS designed for 3d transition metal K edges -> energy range ~5-15 keV
- X-ray tube with silver anode, since Ag characteristic lines above 20 keV

# SBCA - MONOCHROMATING MIRROR

Wavelength selection with crystals  
according to Bragg's law





$$\lambda = 2d_{hkl} \sin \theta_B$$



**Spherically Bent Crystal Analyser**


Monochromatizes *and* focuses the beam

# MOTORIZED JOHANN GEOMETRY






# SCINTILLATOR DETECTOR



- Doped NaI crystal converts X-ray photons to low energy photons, which again are converted to a cascade of electrons
- Electron pulse is collected and transformed into voltage pulse
- Pulses proportional to photon energy but poor resolution (~30%)

# DEAD-TIME CORRECTION




- The detector has a finite response time to incident photons, known as the dead time  $\tau$
- During that time, the detector can not count other photons
- If the effect is fairly small, can be corrected using


$$N_{\text{correct}} \approx \frac{N}{1 - \frac{\tau}{T}N}$$

where  $N$  is the number of counts and  $T$  is the time in which they were acquired

# ABSORPTION IN TRANSMISSION GEOMETRY




# BACKGROUND NOISE



- In addition to signal, there is always background of elastic scattering, fluorescence, etc. present
- $\mu x$  from logarithm = highly non-linear = sensitive to background

# BACKGROUND REMOVAL



- Measure background by moving the detector away from the direct beam
- For accuracy, measure on both sides of the focus and take the mean
- Background slowly varying  $\rightarrow$  fit low-order polynomial and remove the fit from the signal to avoid increased statistical noise



# MEASUREMENT PROCEDURE (ABSORPTION)

1. Measure the direct beam  $I_0$ , the transmitted beam  $I$ , and their backgrounds  $I_{0,\text{bg}}$  and  $I_{\text{bg}}$
2. Apply the dead-time correction
3. Fit the low order polynomials  $y$  and  $y_0$  to the backgrounds  $I_{0,\text{bg}}$  and  $I_{\text{bg}}$
4. Compute  $\mu x$  from the equation:

$$\mu x = \ln \frac{I_0 - y_0}{I - y}$$



# EXAMPLE: Co foil in transmission




# EXAMPLE: Co foil in transmission

Energy scale conversion:

$$\lambda = 2d_{hkl} \sin \theta$$


$$E = \frac{hc}{\lambda}$$






# ABSORPTION THROUGH FLUORESCENCE

Ionization of a core electron  
leaves a hole behind




Hole is filled and a fluorescence  
photon is emitted




$$I_f(E) \propto \mu(E)$$

for thin sample

# FLUORESCENCE MODE



# SOLID STATE DETECTOR



- For fluorescence signal, we need to be able to separate the fluorescence lines -> better energy resolution needed
- CdTe-detector is used with multichannel analyser
- Saturates easily, not good for transmission

# WHEN TO USE FLUORESCENCE?

- Never, if the transmission mode works
- In fluorescence: no absolute absorption coefficient, not good for (optically) thick samples, low photon yield, geometric difficulties due to spherical aberration
- However, useful for dilute samples, thin films on thick substrate, some complex environments



# HelixAS

X-RAY SOURCE

MUCH PHOTON,  
VERY SCIENCE

DETECTORS

ONLY 90 k€\*!

MOTORS AND MOVEMENT

\*Final cost may vary

